



GEYSER CLEANING | DESCRIPTION





One of the most important installation on a dental clinic is the suction system. If the system is not maintained and cleaned every day with a disinfect, fluid blocks could easily occur. Rectus has more than 40 years of experience in production and servicing of dental suction system. Our knowledge throughout all these years have resulted in Geysers - a highly effective fluid specialized for dental suction systems.

Geysers disinfection both de-grease and disinfect the suction system. Using Geysers on a daily basis prevents building up of biofilm and blocks in filters and the pipe system. Furthermore Geysers fluid prevents odor problems.

DIMENSIONS

PRODUCT	ITEM NO.	PACKAGE	WEIGHT
Geysers cleaning	890011	6 pcs.	6 kg

RELATED PRODUCTS

	<table border="1"> <thead> <tr> <th>PRODUCT</th> <th>PAGE NO.</th> </tr> </thead> <tbody> <tr> <td>Separation tanks</td> <td>ST 01</td> </tr> </tbody> </table>	PRODUCT	PAGE NO.	Separation tanks	ST 01		<table border="1"> <thead> <tr> <th>PRODUCT</th> <th>PAGE NO.</th> </tr> </thead> <tbody> <tr> <td>Domega suction motor</td> <td>SM 02</td> </tr> </tbody> </table>	PRODUCT	PAGE NO.	Domega suction motor	SM 02
PRODUCT	PAGE NO.										
Separation tanks	ST 01										
PRODUCT	PAGE NO.										
Domega suction motor	SM 02										
	<table border="1"> <thead> <tr> <th>PRODUCT</th> <th>PAGE NO.</th> </tr> </thead> <tbody> <tr> <td>Sep-it</td> <td>ST 02</td> </tr> </tbody> </table>	PRODUCT	PAGE NO.	Sep-it	ST 02		<table border="1"> <thead> <tr> <th>PRODUCT</th> <th>PAGE NO.</th> </tr> </thead> <tbody> <tr> <td>R5 Suction motor</td> <td>SM 01</td> </tr> </tbody> </table>	PRODUCT	PAGE NO.	R5 Suction motor	SM 01
PRODUCT	PAGE NO.										
Sep-it	ST 02										
PRODUCT	PAGE NO.										
R5 Suction motor	SM 01										

